


Introducing the MADSEN Astera

myaudiometer. It's familiar,
and it's so much


More



MADSEN Astera



otometrics
MADSEN • AURICAL • ICS



**"I'm amazed that
I can do so much
more without
changing the way
I work."**

Janice Audiologist

myaudiometer. It's familiar. Just more...

Possibilities

Clinical audiometry has never been so flexible.

The MADSEN Astera opens up a world of testing possibilities with features and options that let you work the way you want. It's the first clinical audiometer that combines the best of traditional, stand-alone audiometers and the newer PC-based systems. With the Astera, you can choose to operate through the familiar user interface of the Audiometer Control Panel or the PC keyboard or mouse.

The audiometer features a comprehensive test battery that fulfills all your clinical audiometry requirements. To cover the needs of multiple users and various testing scenarios, the Astera can be personalized with the extensive user test functionality, which quickly lets you choose your favorite test settings and sequences. The result is a smoother and faster work process.


Learn more:

www.myaudiometer.com



Audiometer Control Panel

The full function Audiometer Control Panel can handle all possible testing and allows a smooth transition to a PC based system.



Unique user test possibilities

The user test functionality actively supports the clinical workflow and makes testing seamless and efficient.



Integrated speech material

The Astera comes standard with integrated speech material. It significantly improves the testing speed and user comfort while testing.



Customized view options

Astera can easily be adapted to your specific viewing preferences when testing. These preferences can also be implemented in your user tests.



“Finally, a true two-channel clinical audiometer that combines the best of two worlds and delivers high performance all the way from testing to storing and sharing results.”

Shane Research Audiologist

myaudiometer. It's familiar. Just more...

Power

With its extensive test battery and true two-channel setup, the MADSEN Astera fully qualifies as a clinical audiometer.

The proven technical platform meets all ANSI standards and the powerful built-in amplifier ensures high performance during sound-field testing. A unique calibration method that ensures sound accuracy for all frequencies up to 20,000 Hz further strengthens the Astera's capabilities.

The Astera also features a familiar, user-friendly interface and a design that adds a modern twist to true traditional audiometry.

In short, the Astera is simply smarter, sleeker and much easier to work with than anything you've ever seen in clinical audiometry.

Learn more:

www.myaudiometer.com



High frequency testing

The Astera has the ability to perform high frequency testing up to 20,000 Hz – which makes it ideal for tinnitus testing and ototoxic monitoring.



A built-in power amplifier

Astera enables you to perform sound field testing at 90 dB HL in a 6' by 6' sound room.



Advanced patient communication

The Astera comes standard with an advanced patient communication system, including a custom-made monitor headset, TF/TB microphones and a talk-to-assistant functionality.



SNR on the speech screen

Astera allows electronic storage of any Signal to Noise Ratio (SNR) regardless of the source. Also it facilitates the use of contra lateral masking.



True two-channel clinical audiometer

With an extensive clinical test battery, flexible mixing/routing capabilities, and a true two-channel interface, you get full flexibility and a clear overview while testing.





**"It's so nice that the
Audiometer Control
Panel allows me to keep
my eyes on the child
while testing."**

Carrie Pediatric Audiologist



myaudiometer. It's familiar. Just more...

Care

When testing, it is important to be in close contact with the patient. The MADSEN Astera's Audiometer Control Panel features a familiar and intuitive layout that lets you focus on the patient without looking at the instrument. The testing process is furthermore supported by a talk-to-assistant functionality, a children's picture board wordlist, large audiograms and several viewing modes.

In addition to fulfilling all traditional pediatric testing needs, the Astera is also in the lead, when it comes to new testing methods such as using FRESH (FREquency Specific Hearing assessment) noise and a fully flexible speaker routing system.

Learn more:

www.myaudiometer.com



FRESH noise

Otometrics introduces a third stimulus that is appropriately calibrated for pediatric testing; FRESH (FREquency Specific Hearing assessment) noise.



Audiometer Control Panel

An ergonomic, impressive design and an intuitive layout, make the Audiometer Control Panel, very user friendly and efficient to use.



Multiple speaker routing possibilities

You can achieve full flexibility of routing sounds in the testing process. The feature is ideal for behavioral pediatric testing as well as any other aided testing.



Large and clear audiograms

Interpret test data easily and maintain a clear overview while testing.

**“The MADSEN Astera
has increased
our productivity and
improved our workflow
– which makes our staff
and our clients happy.”**

Bob ENT Doctor & Department Director



myaudiometer. It's familiar. Just more...

Efficiency

With the MADSEN Astera, no one in your clinic needs to change their habits – although everyone will enjoy greater freedom, flexibility and the ability to do so much more with just one device.

To meet your specific test requirements and user preferences, you can switch between settings with the press of a button. This ensures easy and quick operation, increases the tester's confidence throughout the process, and reduces testing time.

Since the Astera is part of the OTOsuite universe, the audiometric workflow easily integrates with the immittance testing process. This allows you to print a combined report as well as a wide variety of other professional reports. You can also store data electronically and transfer it to third-party systems – making your clinical workflow even more efficient.


Learn more:
www.myaudiometer.com



Astera is part of the OTOsuite universe
This means a strong integration with other diagnostic and fitting tools – which improves your workflow significantly.




Astera fulfills diverse testing needs
By combining the ability to use the Audiometer Control Panel, a PC keyboard or mouse, the Astera is very efficient when testing – and makes room for individual preferences. This diversity makes the Astera the only audiometer in the industry that can meet all of the audiologist's needs in a busy clinic.



Speed of operation
By using the power of the PC to generate specific user setups and the familiar user interface from the Audiometer Control Panel, the Astera is fast and intuitive to operate.



One-click data logging
By simply clicking a few buttons you can document test findings fast and easy and remove pen and paper from the testing process.



“They say the right tool is half the trick. Given the results we’re seeing, it’s clear we made a smart investment.”

Laura Private Practice Audiologist

myaudiometer. It's familiar. Just more...

Features



- High frequency up to 20,000 Hz
- FRESH (FREquency-Specific Hearing assessment) noise
- Extensive user test functionality
- Automatic PTA calculation for Air and Bone
- Integrated word lists (Q/MASS, picture board list, and VA incl. Maryland CNC, Spanish lists etc.)
- Multiple speaker routing
- Built-in power amplifier
- Masking Assistant™
- Comprehensive hot key access from PC keyboard
- Undo/Redo functionality
- One-click data logging
- Simultaneous two-channel routing
- Combined immittance and audiometry report
- Special tests (ABLB, SISI, Weber, Rinne, Stenger, Tone Decay, Tinnitus tests)
- Talk-to-assistant functionality
- and so much more...




The Audiometer Control Panel allows you full functionality and performance.



The audiometer hardware is small, lightweight and comes standard with a smart cable cover. The dimensions are 12.8 x 10 x 2.4 inches, weight 2.81 lbs.



Descriptive icons on the back of the hardware make it possible to connect headphones, speakers etc. in the right place.



The software has large and clear audiograms and can easily be adjusted to fit your preferred view settings.

myaudiometer. It's familiar. Just more...

Benefits

Fast testing

The full function Astera Audiometer Control Panel ensures high speed of operation throughout the testing process.

Familiarity and comfort

The ergonomic and familiar design of the Audiometer Control Panel provides a test setup that you know and trust and makes the testing process extremely comfortable.

Save your favorite user tests

With the Astera, setting up for tests is faster than ever since you can create, save and recall your preferences for the ideal testing scenario at the push of a button.

Ideal for pediatric testing

Features such as multiple speaker routing, FRESH (FREquency

Specific Hearing assessment) noise and talk-to-assistant functionality, support the pediatric testing process.

Exceptional amplification

The built-in internal power amplifier enables you to perform sound field testing at 90dB without an external amplifier – saving your clinic both money and space.

Total routing flexibility

The true two-channel layout of the software and Audiometer Control Panel give testers full routing flexibility making test setup simple and efficient.

The power to do it all

The Astera's full featured clinical test battery gives you the freedom to perform every type of test necessary in a clinical setting.

Uninterrupted communication

A custom-made monitor headset and Talk-forward/Talk-back microphones, and a talk-to-assistant communication system, ensure clear communication throughout the testing process.

Speech Testing – your way

Integrated word lists in the speech testing module ensure more efficient testing. With the Astera, you see the words on the PC monitor and can set up the speech testing the way you like, so you can, for example, choose to present only certain words.

Focus on what matters

One-click data logging lets you eliminate pen and paper from the testing process – this means a more efficient workflow, better

data handling and more time with your patients.

Personalize your view settings

The Astera allows you to choose your favorite view settings so you easily maintain an overview while testing.

Save your data anywhere

With the Astera it's possible to save data directly on your PC, in NOAH, or in a third party EMR system.

Fast and easy reporting

When you finish testing, you can choose to print the data on a wide variety of different reports, including a combined immittance and audiometry report. This makes documentation fast and easy – and enables you to easily share reports with others.

Small footprint

The audiometer hardware can easily be mounted on the wall, therefore it doesn't take up much space in the clinic.

Enhanced workflow with OTOSuite

Because the Astera audiometer is an integrated part of the OTOSuite universe, it gives you many opportunities to make your workflow more efficient – and even go paperless. OTOSuite combines individual diagnostic and fitting tools in one software and enables you to move seamlessly from one audiometric process to another.

Read more about OTOSuite on www.myaudiometer.com

myaudiometer. It's familiar. Just more...

Value

The MADSEN Astera is directly integrated with the fitting process in the OTOsuite universe. This allows you to perform audiometry, fitting, counseling, and verification in a seamless workflow. You can also store, share, and report all data, as well as connect to third-party systems such as NOAH and EMR.


With smart tools for data registration, integrated word lists and more, the Astera makes pen and paper unnecessary throughout the testing process - so your office can

go paperless! It also enables you to choose between different test reports that can, for example, be personalized with your logo.

As the new leader of true clinical audiometry, the Astera gives you excellent value for your money, improves your professional appearance, and takes your work processes to the next level.

Learn more:

www.myaudiometer.com




Audiogram overlays
With a wide variety of audiogram overlays the Astera enhances the possibilities for patient counseling.




Improved workflow
The clinic's workflow is improved as Astera connects to both counseling and fitting tools in the OTOsuite universe.



Professional reports
Astera offers a wide variety of professional reports, including a combined immittance and audiogram report.



Connection to third party systems
Data can easily be transferred to third party systems such as NOAH and/or Electronic Medical Records.



SII, PTA, masking assistant etc.
A number of specific functionalities make the testing process fast and convenient, as well as ensuring valid test results.

myaudiometer. It's familiar. Just...

Try it

To learn more about the MADSEN Astera
or to book a demonstration go to
www.myaudiometer.com

- **Learn more about the Astera**
- **Take a tour of the Astera**
- **Book a demonstration**
- **The Astera story**

● Hearing Assessment ● Fitting Systems ● Balance Assessment

GN Otometrics,
Europe +45 45 75 55 55,
info@gnotometrics.dk

GN Otometrics,
North America +800 289 2150,
sales@gnotometrics.com

www.otometrics.com



otometrics

MADSEN • AURICAL • ICS